# INSTRUCTIONS FOR USE

## Walker

Product code: **Easy-Walker**

## Short walker

Product code: **Easy-Walker Low**

### 1. Characteristics

The walker consists of a splint in the shape of a shoe which is made of PVC, and inner liner of soft material. Owing to the application of modern materials, the orthosis is air permeable and removes excessive heat and moisture perfectly. The walker is equipped with system of circumferential straps which enable individual adjustment of stabilization force and with profiled anti-sliding sole which facilitates the patient possibly physiological gait. The height of the walker is 17" (Easy-Walker) or 11" (Easy-Walker Low).

### 2. Application

Walker is used in post-traumatic and post-operative rehabilitation process in order to limit the range of mobility in ankle joint and to improve functional activities. The orthosis stabilizes the foot and shin in intermediate (neutral) position, in post-traumatic injuries in the area of ankle joint.

### Indications for use:

- Post-operative treatment
- Twist and dislocation of ankle joint
- Achilles tendon rupture
- As an alternative for plaster cast immobilization
- Fractures in the area of foot, upper and lower ankle joint

### 3. Contraindications

Do not use in case of:

- Skin abrasion
- Eczema
- Wound and other skin injuries
- After application of warming gel or ointment

**Note!** In case of varicose veins the usage of the walker is restricted. Contact with attending physician to determine whether the product may be used together with compression stockings.

### 4. Notes and recommendations

**Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.**

### 5. Fitting instruction

(in sitting position with leg slightly bent in knee joint)

- Unfasten the circumferential straps of the orthosis (velcro straps)
- Pull out the inner liner and unfasten its velcro straps
- Put the inner liner on the leg, and if necessary, use the additional filling cushions included in the set
- Fasten the velcro of the inner liner, then put the inner liner into the orthosis
- Fasten the circumferential straps of the orthosis (velcro straps)
- Check the correctness of the fastened straps in the standing position, and if necessary adjust the strength of their fastening.

#### Size data (shoe size)

S	36-39
M	40-43
L	44-46,5

## 6. Maintenance

- Circumferential straps and inner liner can be washed by hand or in washing machine at temperature of 30°C. Do not tumble dry, chlorinate, iron, bleach, dry clean.
- Ensure the velcro is fastened in order not to destroy the material.
- Dry at room temperature far from intensive heat source.
- Clean the PVC splint with damp cloth with soap.
- Dry the splint thoroughly.

## 7. Storage and transport conditions

Device should not be exposed to direct sunlight, water or moisture.

## 8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

## 9. Labels

	Reference number		Manufacturer
	LOT number		Manufacture date
	Serial number		Medical device
	Note		Please read the instruction
	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: [timago.com](http://timago.com). Thank you for choosing Timago!

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